


# CONSERVATION OF MOMENTUM IN CASE OF EQUAL MASSES

MED 08.06



## Material

| Item-no. | Qty. | Description                      |
|----------|------|----------------------------------|
| DS101-3B | 1    | Stand rail with scale, L=1000 mm |
| DM300-2A | 2    | Dynamics trolley, demo, 50 g     |
| P1312-2A | 2    | Car body for trolley SE          |
| P1311-2D | 1    | Spring bumper                    |
| P3120-2Z | 1    | Universal timer "inno"           |
| P3120-5B | 1    | S-shaped assembly platform       |
| P1320-4A | 2    | Light gate "demo" 04             |
| P1321-3K | 2    | Block for light gates            |
| DS103-1H | 1    | Holder for guide rail            |
| DM362-1E | 1    | Baffle block                     |
| DM344-1S | 1    | Projectile launcher 02           |

# CONSERVATION OF MOMENTUM IN CASE OF EQUAL MASSES

MED 08.06

## Purpose

To demonstrate the conservation of momentum

## Preparation

Position the two blocks for light gates at a distance of 90 cm on the table and place the light gates in between the blocks as shown on the image to the left.

Now place the stand rail on the blocks for light gate.


Mount the projectile launcher on the left end of the stand rail. Afterwards mount the holder for guide rail on the right end of the stand rail and fix the baffle block on the holder.



Place the car bodies on the dynamics trolleys; make sure that the small rod of the trolley points through the hole of the car body.



Attach the spring bumper to the left trolley.



Initially the launching piston of the projectile launcher is fully extended; place the dynamics trolley at the end of the launching piston.


Position the left light gate so that the diodes of it are just in front of the dynamics trolley (approximately at the 32 cm mark).

The right light gate is placed at the 65 cm mark.

Push the launching button of the projectile launcher in and pull up the launching piston, fix the piston at the 4<sup>th</sup> notch.

Place the universal timer on the S-shaped assembly platform.

Connect the light gate with the "L1"-socket of the universal timer and set the switch on the universal timer to "L1 Gate"; afterwards turn the timer on.




Move your finger through the light gate to check if it is working properly.

When passing through with the finger the red LED has to flash up – this time of darkening is measured by the timer. In case that this is not working properly the luminosity has to be adjusted (check at the end of this experiment how this can be done).

# CONSERVATION OF MOMENTUM IN CASE OF EQUAL MASSES

MED 08.06

## Experiment




Push the „Reset“ -button on the timer; afterwards push the dynamics trolley completely against the launching piston. By pressing the launching button the dynamics trolley gets pushed away.

The universal timer measures the "darkening time" - the time in which the trolley moves through the laser barrier of the light gate.



Based on the length of the dynamics trolley (125 mm) and the darkening time the current speed of the dynamics trolley can be calculated:

$$0.125 \text{ m} / \dots\dots\dots \text{ s} = \dots\dots\dots \text{ m/s}$$



Again the launching piston is pulled up and fixed at the 4th notch as before; the dynamics trolleys are placed in their original position again.

Instead of the left light gate connect the right light gate now to the "L1"-socket of the timer.

# CONSERVATION OF MOMENTUM IN CASE OF EQUAL MASSES

MED 08.06



Push the „Reset“ -button on the timer again.  
Push the dynamics trolley completely against the launching piston. By pressing the launching button the dynamics trolley gets pushed away.

Now the universal timer measures the "darkening time" at the right light gate.


As before the current speed of the dynamics trolley is calculated again:

$$0.125 \text{ m} / \dots\dots\dots \text{ s} = \dots\dots\dots \text{ m/s}$$

### Result

The pushing trolley comes fully to rest, its momentum becomes Zero.  
The pushed trolley experiences the same force (in opposite direction) and takes over the energy and the momentum since it is accelerated due to the equal masses to the same speed that the other trolley had before it got braked.

If the pushing trolley has a masse  $m$  and the speed  $v$  its momentum is  $p = m \times v$ .



### Note

We will find out that the speed of the second trolley is slightly lower, this loss is created by the friction of the trolley and through conversing the energy at the impact.

### Adjusting the sensitivity of the light gate

Turn the regulator clockwise until it does not go any further; the LED must light up now.  
Afterwards the regulator is slowly turned counterclockwise until the LED no longer lights up.  
The sensitivity is at the highest in this setting.  
The adjustment must be carried out especially when experimenting with very high or very low room brightness.

# CONSERVATION OF MOMENTUM IN CASE OF EQUAL MASSES

MED 08.06

Quantitative evaluations

The trolleys form a closed system, thus the principle of energy and linear momentum applies.

The following equations for the speeds can be derived:

$$v_1' = \frac{v_1 (m_1 - m_2) + 2m_2 \times v_2}{m_1 + m_2}$$

$$v_2' = \frac{v_2 (m_2 - m_1) + 2m_1 \times v_1}{m_1 + m_2}$$

In our experiment the masses are equal, therefore:

$$v_2' = \frac{0 + 2m_1 \times v_1}{2m_1} = v_1$$

$$v_1' = \frac{0 + 2m_2 \times 0}{2m_2} = 0$$