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
The goniometer consists of a base with a rail mounted on it, as seen in the image on the right. In the middle of the base is a pivot axis, on which the left rail and angular scale are moved at the middle.

Function of the screws: the screw on the left side (**S1**) releases the movable arm with the rail, while the right front screw with the smaller diameter (**S2**) serves to fix the angular scale. The "big" front screw on the left (**S3**) allows for fine adjustment of the rail, but only when the other two screws have been tightened. Since the receiver does not have a narrow aperture, setting to a fraction of a degree is not necessary during ultrasound testing; i.e., the screw does not need to be used.

Basic setting:

1. Loosen the screws for the movable arm S1 and the angular scale S2.
2. Align the two rails with a long ruler, a long rail ... or "by sight". Then, secure it with screw S1.
3. Set the angular scale to 0° (the pointer of the vernier scale, which is connected to the arm, points to 0°).
4. Secure the angular scale with the right screw S2.

Now, when the left rail is moved, the arrow on the vernier scale shows the rotation angle (or we have to calculate the difference to the angle after alignment).

With the rail in the vertical position on the angular scale:

Set the pointer to 90° as in no. 3 above.